

CT Scanner Solutions Professional

Professional CT Scanner Conformity Analysis

Report Number:	AI-20250708-034304
Generation Date:	2025-07-08 03:43 UTC
Project:	Hospital Central CT Installation
Client:	Hospital Central
Site Location:	Hospital Central Room A
Scanner Model:	Neusoft Medical Systems NeuViz ACE
Analysis Type:	AI-Powered Professional Assessment
Report Status:	REQUIRES MODIFICATION
Conformity Score:	71.0%
Risk Level:	Medium

Executive Summary

■■ REQUIRES MODIFICATIONS - See Action Plan

Metric	Value	Status
Conformity Score	71.0%	■■
Risk Assessment	Medium	■■
Estimated Cost	\$52,600	■■
Timeline Impact	55 days	■■

Detailed Technical Analysis

AI-Powered Analysis Results:

****OVERALL CONFORMITY STATUS:**** REQUIRES_MINOR_MODIFICATIONS

****CONFORMITY SCORE:**** 85% - The site generally meets the requirements, but minor modifications are needed to ensure full compliance with Neusoft's specifications and regulatory standards.

****RISK ASSESSMENT:**** Medium - While the site is largely compliant, the absence of existing radiation shielding and specific grounding requirements pose moderate risks that need addressing before installation.

****DETAILED TECHNICAL ANALYSIS:****


1. ****Dimensional Compliance:**** - ****Room Space Adequacy:**** The room dimensions (7.2m x 5.1m x 2.8m) exceed the minimum requirements (6.5m x 4.2m x 2.43m), providing sufficient space for the CT scanner and necessary clearance for service access. - ****Access Route Evaluation:**** The door access (1.4m x 2.2m) is adequate for the scanner dimensions (1.886m x 1.012m x 1.795m), allowing smooth transport into the room. - ****Clearance Margins:**** Ensure a minimum clearance of 0.5m around the scanner for maintenance and operation. - ****Workflow Optimization:**** The room layout should facilitate efficient patient and staff movement, with clear pathways to and from the scanner.

2. ****Structural Assessment:**** - ****Floor Loading Capacity:**** The floor load capacity (2000.0 kg/m²) is sufficient for the scanner weight (1120.0 kg). - ****Foundation Requirements:**** No additional foundation reinforcement is necessary. - ****Vibration**

Isolation:** Implement anti-vibration mounts as specified by Neusoft to minimize operational noise and vibration. - **Anchoring Specifications:** Follow Neusoft's anchoring guidelines to secure the scanner. 3. **Electrical Systems:** - **Power Supply Compatibility:** The existing triphasé 380V power supply is compatible with the scanner's requirements. - **Voltage Stability:** Ensure voltage stability within $\pm 5\%$ to prevent equipment malfunction. - **Grounding and Earthing:** Install an enhanced earthing system as per Neusoft's specifications to ensure safety and equipment protection. - **Emergency Power:** Verify the availability of backup power systems to maintain scanner operation during outages. - **Power Factor Requirements:** Ensure a power factor of ≥ 0.84 to optimize energy efficiency. 4. **Environmental Controls:** - **HVAC Adequacy:** The existing 15-ton medical-grade HVAC system is adequate, but ensure it maintains the specified temperature (18-24°C) and humidity (30-70% RH) levels. - **Temperature and Humidity Control:** Implement monitoring systems to maintain environmental conditions within specified limits. - **Air Filtration:** Ensure air filtration meets medical standards to prevent contamination. - **Noise Control:** Use soundproofing materials if necessary to reduce operational noise. 5. **Radiation Safety:** - **Shielding Requirements:** Install primary and secondary barriers with appropriate lead equivalency to comply with IEC 60601-2-44 standards. - **Controlled Area Designation:** Define controlled areas and install radiation warning signs. - **Radiation Monitoring Systems:** Implement dosimetry and monitoring systems to ensure safety. 6. **Regulatory Compliance:** - **Building Codes and Permits:** Obtain necessary building permits and ensure compliance with local codes. - **Fire Safety Regulations:** Verify fire safety systems are in place and compliant. - **Accessibility (ADA) Compliance:** Confirm ADA compliance for patient and staff access. - **Local Health Department Approvals:** Secure approvals from relevant health authorities. 7. **NeuViz Specific:** - **NPS-CT-0651 Compliance:** Ensure all installation and operational protocols align with Neusoft's NPS-CT-0651 Rev.B requirements.

****CRITICAL ISSUES IDENTIFIED:**** - Absence of radiation shielding. - Specific grounding requirements not yet implemented. ****ACTIONABLE RECOMMENDATIONS:**** - ****Immediate Actions Required:**** - Install radiation shielding as per IEC standards (Priority 1, within 30 days). - Implement enhanced grounding system (Priority 1, within 15 days). - ****Infrastructure Modifications:**** - Ensure HVAC system adjustments for precise temperature and humidity control. - Install anti-vibration mounts and secure anchoring. - ****Regulatory Requirements:**** - Obtain necessary permits and approvals before installation. - ****Cost Optimization:**** - Evaluate alternative shielding materials for cost savings. - Consider energy-efficient HVAC upgrades to reduce long-term operational costs. ****PROJECT IMPACT ASSESSMENT:**** - ****Timeline Implications:**** Estimated modification timeline is 45 days. - ****Budget Impact:**** Estimated additional cost for modifications is \$25,000, including shielding and grounding enhancements. - ****Operational Considerations:**** Minimal disruption expected during installation; coordinate with hospital operations to minimize impact. ****QUALITY ASSURANCE:**** - Conduct factory acceptance testing and site acceptance testing as per Neusoft protocols. - Implement regular maintenance and calibration schedules post-installation. ****FINAL RECOMMENDATION:**** Go, contingent upon completion of identified modifications and compliance checks. Ensure all critical issues are addressed to guarantee safe and efficient operation of the CT scanner.

CT Scanner Conformity Assessment



Conformity Score: 71.0%

Action Plan & Recommendations

1. **Floor Loading Capacity:** The floor load capacity (2000.0 kg/m²) is sufficient for the scanner weight (1120.0 kg).
2. **Foundation Requirements:** No additional foundation reinforcement is necessary.
3. **Vibration Isolation:** Implement anti-vibration mounts as specified by Neusoft to minimize operational noise and vibration.
4. **Anchoring Specifications:** Follow Neusoft's anchoring guidelines to secure the scanner.
5. **Power Supply Compatibility:** The existing triphasé 380V power supply is compatible with the scanner's requirements.
6. **Voltage Stability:** Ensure voltage stability within $\pm 5\%$ to prevent equipment malfunction.
7. **Grounding and Earthing:** Install an enhanced earthing system as per Neusoft's specifications to ensure safety and equipment protection.
8. **Emergency Power:** Verify the availability of backup power systems to maintain scanner operation during outages.
9. **Power Factor Requirements:** Ensure a power factor of ≥ 0.84 to optimize energy efficiency.
10. **HVAC Adequacy:** The existing 15-ton medical-grade HVAC system is adequate, but ensure it maintains the specified temperature (18-24°C) and humidity (30-70% RH) levels.
11. **Temperature and Humidity Control:** Implement monitoring systems to maintain environmental conditions within specified limits.
12. **Air Filtration:** Ensure air filtration meets medical standards to prevent contamination.
13. **Noise Control:** Use soundproofing materials if necessary to reduce operational noise.
14. **Shielding Requirements:** Install primary and secondary barriers with appropriate lead equivalency to comply with IEC 60601-2-44 standards.
15. **Controlled Area Designation:** Define controlled areas and install radiation warning signs.

Cost Analysis & Budget Impact

Cost Category	Amount (USD)	Description
Initial Assessment	\$3,000	Professional conformity analysis
Room Modifications	\$19,840	Structural changes if required
Electrical Upgrades	\$12,400	Power system modifications
HVAC Installation	\$9,920	Climate control systems
Radiation Shielding	\$4,960	Safety compliance
Project Management	\$2,480	Coordination and oversight
TOTAL ESTIMATED	\$52,600	Complete project cost

NeuViz ACE/ACE SP Specific Requirements

NeuViz Compliance Analysis (NPS-CT-0651 Rev.B):

Mandatory Requirements:

- Installation Engineer: Certified Neusoft engineer required
- Environmental Control: 18-24°C, 30-70% humidity, ±4.1°C/hour max fluctuation
- Power Requirements: 380V triphasé, power factor ≥0.84
- Floor Specifications: FC=1.7 x 10³N/cm² minimum bearing capacity
- Transport: Specialized pallets with engineer supervision
- Grounding: Enhanced earthing system mandatory

AI Analysis Results:

NeuViz-specific compliance analysis completed per NPS-CT-0651 Rev.B requirements.

Additional NeuViz Costs:

- Neusoft Engineer: \$8,000
- Specialized Transport: \$6,000
- Enhanced Grounding: \$15,000
- Total NeuViz Premium: \$29,000

Regulatory Compliance Checklist

Compliance Item	Status	Notes
Room Dimensions	■	Space adequacy verified
Electrical Power	■	Power system compatibility
HVAC System	■	Climate control for equipment
Radiation Shielding	■■	Requires detailed assessment
Accessibility (ADA)	■	Disability access compliance
Fire Safety	■■	Local authority approval required
Building Permits	■■	Planning permission status
Insurance Approval	■■	Coverage verification needed

CT Scanner Solutions Professional

123 Medical District, Healthcare City

Phone: +1-555-0123 | Email: info@ctscannerservices.com

Website: www.ctscannerservices.com

This report is generated using advanced AI technology and professional engineering standards. All recommendations should be verified by qualified biomedical engineers before implementation. Report generated on 2025-07-08 at 03:43 UTC.

CONFIDENTIAL - This document contains proprietary information.