CT Scanner Solutions Professional Professional CT Scanner Conformity Analysis

Report Number:	AI-20250708-034316
Generation Date:	2025-07-08 03:43 UTC
Project:	Hospital Central CT Installation
Client:	Hospital Central
Site Location:	Hospital Central Room A
Scanner Model:	Neusoft Medical Systems NeuViz ACE
Analysis Type:	AI-Powered Professional Assessment
Report Status:	REQUIRES MODIFICATION
Conformity Score:	71.0%
Risk Level:	Medium

Executive Summary

REQUIRES MODIFICATIONS - See Action Plan

Metric	Value	Status
Conformity Score	71.0%	
Risk Assessment	Medium	
Estimated Cost	\$52,600	
Timeline Impact	55 days	


Detailed Technical Analysis

Al-Powered Analysis Results:

OVERALL CONFORMITY STATUS: REQUIRES MINOR MODIFICATIONS **CONFORMITY SCORE:** 85% - The installation site meets most of the requirements, but minor modifications are needed for full compliance, particularly in electrical grounding and radiation shielding. **RISK ASSESSMENT:** Medium - The primary risks involve electrical grounding adequacy and radiation shielding, which can be addressed with targeted interventions. **DETAILED TECHNICAL ANALYSIS:** 1. **Dimensional Compliance:** - **Room Space Adequacy:** The room dimensions (7.2m x 5.1m x 2.8m) exceed the minimum requirements (6.5m × 4.2m × 2.43m), providing sufficient space for the CT scanner and necessary clearance for service access. - **Access Route Evaluation:** The door access (1.4m × 2.2m) is adequate for the scanner dimensions (1.886 x 1.012 x 1.795), allowing smooth transport and installation. - **Clearance Margins:** Adequate clearance is available for maintenance and service access, ensuring optimal workflow for patients and staff. 2. **Structural Assessment:** - **Floor Loading Capacity:** The floor load capacity (2000.0 kg/m²) is sufficient for the scanner weight (1120.0 kg), meeting the minimum bearing capacity requirement of FC=1.7 x 103N/cm². **Vibration Isolation:** Specialized anti-vibration mounting and anchoring systems are recommended to minimize operational vibrations. 3. **Electrical Systems:** - **Power

Supply Compatibility:** The existing triphasé 380V power supply is compatible with the scanner's requirements. - **Grounding and Earthing:** Enhanced earthing systems are required to meet specialized grounding specifications. This may involve additional grounding rods or plates. - **Emergency Power:** Ensure integration with hospital's emergency power systems to maintain operation during outages. - **Power Factor:** Ensure power factor correction to meet the ≥0.84 requirement. 4. **Environmental Controls:** - **HVAC Adequacy:** The existing 15-ton medical-grade HVAC system is adequate for cooling requirements. Ensure air conditioning does not directly impact the patient table. - **Temperature and Humidity Control:** The system must maintain 18-24°C and 30-70% RH, with fluctuations not exceeding ±4.1°C/hour. 5. **Radiation Safety:** -**Shielding Requirements:** As existing shielding is absent, install primary and secondary barriers with appropriate lead equivalency based on IEC 60601-2-44 standards. -**Controlled Area Designation:** Define and mark controlled areas to ensure radiation safety compliance. 6. **Regulatory Compliance:** - **Permits and Codes:** Verify compliance with local building codes and obtain necessary permits. - **Fire Safety and Accessibility:** Ensure fire safety measures and ADA compliance are maintained. 7. **NeuViz Specific:** - **NPS-CT-0651 Compliance:** Ensure all environmental and power requirements are met, including specialized transport and installation supervision by a certified Neusoft engineer. **CRITICAL ISSUES IDENTIFIED:** - Inadequate grounding system. - Absence of radiation shielding. **ACTIONABLE RECOMMENDATIONS:** -**Immediate Actions Required:** - Upgrade grounding system to meet enhanced requirements. - Install radiation shielding as per IEC standards. - **Infrastructure Modifications:** - Implement anti-vibration mounts and ensure proper air circulation. -**Regulatory Requirements:** - Obtain necessary permits and approvals from local health authorities. - **Cost Optimization:** - Evaluate alternative grounding solutions to reduce costs. **PROJECT IMPACT ASSESSMENT:** - **Timeline Implications:** Estimated modification timeline is 30-45 days, considering procurement and installation of shielding and grounding systems. - **Budget Impact:** Estimated additional costs for grounding and shielding are \$15,000-\$20,000. - **Operational Considerations:** Minimal disruption expected during installation; coordinate with hospital operations to minimize impact. **QUALITY ASSURANCE:** - Conduct factory acceptance testing and site acceptance testing as per Neusoft protocols. **FINAL RECOMMENDATION:** Go - Proceed with installation after addressing identified modifications, ensuring all compliance and safety standards are met.

CT Scanner Conformity Assessment



Conformity Score: 71.0%

Action Plan & Recommendations

- 1. **Emergency Power:** Ensure integration with hospital's emergency power systems to maintain operation during outages.
- 2. **Power Factor:** Ensure power factor correction to meet the ≥0.84 requirement.
- 3. **HVAC Adequacy:** The existing 15-ton medical-grade HVAC system is adequate for cooling requirements. Ensure air conditioning does not directly impact the patient table.
- 4. **Shielding Requirements:** As existing shielding is absent, install primary and secondary barriers with appropriate lead equivalency based on IEC 60601-2-44 standards.
- 5. **Controlled Area Designation:** Define and mark controlled areas to ensure radiation safety compliance.
- 6. **Permits and Codes:** Verify compliance with local building codes and obtain necessary permits.
- 7. **Fire Safety and Accessibility:** Ensure fire safety measures and ADA compliance are maintained.
- 8. **NPS-CT-0651 Compliance:** Ensure all environmental and power requirements are met, including specialized transport and installation supervision by a certified Neusoft engineer.

Cost Analysis & Budget Impact

Cost Category	Amount (USD)	Description
Initial Assessment	\$3,000	Professional conformity analysis
Room Modifications	\$19,840	Structural changes if required
Electrical Upgrades	\$12,400	Power system modifications
HVAC Installation	\$9,920	Climate control systems
Radiation Shielding	\$4,960	Safety compliance

Project Management	\$2,480	Coordination and oversight
TOTAL ESTIMATED	\$52,600	Complete project cost

NeuViz ACF/ACF SP Specific Requirements

NeuViz Compliance Analysis (NPS-CT-0651 Rev.B):

Mandatory Requirements:

- Installation Engineer: Certified Neusoft engineer required
- Environmental Control: 18-24°C, 30-70% humidity, ±4.1°C/hour max fluctuation
- Power Requirements: 380V triphasé, power factor ≥0.84
- Floor Specifications: FC=1.7 x 10³N/cm² minimum bearing capacity
- Transport: Specialized pallets with engineer supervision
- Grounding: Enhanced earthing system mandatory

Al Analysis Results:

NeuViz-specific compliance analysis completed per NPS-CT-0651 Rev.B requirements.

Additional NeuViz Costs:

Neusoft Engineer: \$8,000
Specialized Transport: \$6,000
Enhanced Grounding: \$15,000
Total NeuViz Premium: \$29,000

Regulatory Compliance Checklist

Compliance Item	Status	Notes
Room Dimensions		Space adequacy verified
Electrical Power		Power system compatibility
HVAC System		Climate control for equipment
Radiation Shielding		Requires detailed assessment
Accessibility (ADA)		Disability access compliance
Fire Safety		Local authority approval required
Building Permits		Planning permission status
Insurance Approval		Coverage verification needed

CT Scanner Solutions Professional

123 Medical District, Healthcare City

Phone: +1-555-0123 | Email: info@ctscannerservices.com

Website: www.ctscannerservices.com

This report is generated using advanced AI technology and professional engineering standards. All recommendations should be verified by qualified biomedical engineers before implementation. Report generated on 2025-07-08 at 03:43 UTC.