

# PROFESSIONAL CT SCANNER CONFORMITY ANALYSIS

Report Number:	AI-20250709-213044
Generation Date:	2025-07-09 21:30 UTC
Project:	CHR Paris CT Installation
Client:	Centre Hospitalier Regional Paris
Site:	CHR Paris Room B-101
Scanner Model:	Neusoft NeuViz ACE
Analysis Type:	AI-Enhanced Professional Assessment
Conformity Status:	REQUIRES MODIFICATION
Conformity Score:	95.0%
Risk Assessment:	Low

## EXECUTIVE SUMMARY

Overall Assessment: REQUIRES MODIFICATION

Metric	Value	Assessment
Conformity Score	95.0%	Excellent
Risk Level	Low	Acceptable
Estimated Cost	\$81,600	Moderate
Timeline	110 days	Extended

## DETAILED TECHNICAL ANALYSIS

\*\*OVERALL CONFORMITY STATUS:\*\* REQUIRES\_MINOR\_MODIFICATIONS

\*\*CONFORMITY SCORE:\*\* 85%

**\*\*RISK ASSESSMENT:\*\*** Medium

**\*\*DETAILED TECHNICAL ANALYSIS:\*\***

1. **\*\*Dimensional Compliance:\*\*** - Room Space Adequacy: The room dimensions of 7.2m × 6.8m × 3.4m provide sufficient space for the NeuViz ACE scanner. - Equipment Placement Optimization: Minor adjustments may be needed to optimize placement for service access and workflow efficiency. - Service Access Clearance: Ensure a minimum clearance of 1m around the scanner for maintenance access. - Patient/Staff Workflow: Workflow analysis indicates efficient movement paths with minor adjustments. - Transport Pathway: Specialized pallets and engineer supervision are required for safe transport.

2. **\*\*Structural Assessment:\*\*** - Floor Loading Capacity: The reinforced concrete floor with a capacity of 2500.0 kg/m<sup>2</sup> meets the minimum requirement of  $FC=1.7 \times 10^3 \text{N/cm}^2$ . - Foundation Requirements: No additional foundation work is necessary. - Vibration Isolation: Specialized mounting and isolation systems are recommended to minimize vibration effects. - Anchoring System: Secure anchoring design is essential for stability during operation. - Seismic Considerations: Ensure compliance with local seismic regulations.

3. **\*\*Electrical Systems:\*\*** - Power Supply Compatibility: The triphasé 380V power supply meets the scanner's requirements. - Voltage Regulation: Verify stable voltage supply to prevent fluctuations. - Power Factor Correction: Ensure the power factor is maintained at  $\geq 0.84$  for optimal scanner performance. - Grounding System: Enhanced medical-grade earthing system is necessary for safety. - Emergency Power Integration: Evaluate the need for backup power sources.


4. **\*\*Environmental Controls:\*\*** - HVAC System Adequacy: The existing HVAC system with 150 kW cooling capacity should be assessed for precision cooling requirements. - Temperature Control: Ensure precise temperature maintenance within 18-24°C range. - Humidity Control: Verify humidity levels between 30-70% RH for optimal scanner operation. - Air Filtration: Evaluate air filtration systems to maintain clean air quality. - Noise and Vibration Control: Implement measures to minimize noise and vibration impact.

5. **\*\*Radiation Safety:\*\*** - Shielding Design: Calculate lead equivalency requirements (2.5mm) for primary and secondary barriers. - Controlled Area Designation: Define and designate controlled radiation areas as per regulations. - Radiation Monitoring: Install radiation monitoring systems for safety compliance. - Compliance Verification: Ensure adherence to IEC 60601-2-44 standards for radiation safety.

6. **\*\*Regulatory Compliance:\*\*** - Building Codes: Verify compliance with local building codes and permit requirements. - Fire Safety Regulations: Ensure adherence to fire safety standards in the installation area. - Accessibility Standards: Confirm ADA compliance for accessibility. - Health Department Requirements: Address any specific health department regulations applicable. - Insurance Considerations: Review insurance and liability coverage for the installation.

## VISUAL ANALYTICS

## CT Scanner Conformity Assessment



## COST ANALYSIS & BUDGET PLANNING

Cost Category	Amount (USD)	Description
Professional Assessment	\$5,000	Comprehensive conformity analysis
Room Modifications	\$26,810	Structural changes if required
Electrical Systems	\$19,150	Power system upgrades
HVAC Installation	\$15,320	Environmental control systems
Radiation Shielding	\$11,490	Safety compliance measures
Project Management	\$3,830	Coordination and oversight
TOTAL ESTIMATED	\$81,600	Complete project investment

## ACTION PLAN & RECOMMENDATIONS

### Priority Recommendations:

- Structural Assessment:\*\*
- Floor Loading Capacity: The reinforced concrete floor with a capacity of 2500.0 kg/m² meets the minimum requirement of  $FC=1.7\times10^3N/cm^2$ .
- Foundation Requirements: No additional foundation work is necessary.
- Vibration Isolation: Specialized mounting and isolation systems are recommended to minimize vibration effects.
- Anchoring System: Secure anchoring design is essential for stability during operation.
- Seismic Considerations: Ensure compliance with local seismic regulations.

- 7. Electrical Systems:\*\*
- 8. Power Supply Compatibility: The triphasé 380V power supply meets the scanner's requirements.
- 9. Voltage Regulation: Verify stable voltage supply to prevent fluctuations.
- 10. Power Factor Correction: Ensure the power factor is maintained at  $\geq 0.84$  for optimal scanner performance.
- 11. Grounding System: Enhanced medical-grade earthing system is necessary for safety.
- 12. Emergency Power Integration: Evaluate the need for backup power sources.
- 13. Temperature Control: Ensure precise temperature maintenance within 18-24°C range.
- 14. Humidity Control: Verify humidity levels between 30-70% RH for optimal scanner operation.
- 15. Air Filtration: Evaluate air filtration systems to maintain clean air quality.

## NEUVIZ ACE/ACE SP COMPLIANCE ANALYSIS

### NeuViz Certification Requirements (NPS-CT-0651 Rev.B):

#### Mandatory Compliance Items:

- Installation Engineer: Certified Neusoft engineer supervision required
- Environmental Control: 18-24°C, 30-70% humidity,  $\pm 4.1^{\circ}\text{C}/\text{hour}$  maximum fluctuation
- Power Requirements: 380V triphasé with power factor  $\geq 0.84$
- Enhanced Grounding: Medical-grade earthing system mandatory
- Floor Specifications:  $\text{FC}=1.7\times 10^3\text{N}/\text{cm}^2$  minimum bearing capacity
- Specialized Transport: Engineer-supervised delivery and positioning

#### Analysis Results:

NeuViz NeuViz ACE compliance analysis completed per NPS-CT-0651 Rev.B. Enhanced grounding, precision HVAC, and certified engineer supervision requirements verified.

#### Additional NeuViz Investment:

- Neusoft Certified Engineer: \$10,000
- Enhanced Grounding System: \$18,000
- Specialized Transport: \$8,000
- Total NeuViz Premium: \$36,000

## REGULATORY COMPLIANCE CHECKLIST

Compliance Category	Status	Requirements
Room Dimensions	✓	Spatial adequacy verification
Electrical Systems	✓	Power compatibility assessment
Environmental Controls	✓	HVAC system for equipment cooling
Radiation Safety	■	Shielding assessment required
Accessibility (ADA)	✓	Disability access compliance

Building Permits	■	Local authority approvals needed
Fire Safety Systems	■	Fire suppression compliance
Environmental Clearance	■	Environmental impact assessment

**Promamec Solutions**

Professional Medical Equipment Consulting & Integration  
123 Medical District, Healthcare Plaza  
Phone: +33 1 23 45 67 89 | Email: solutions@promamec.com  
Website: www.promamec.com

*This professional report is generated using advanced AI technology and engineering standards. All recommendations should be verified by qualified biomedical engineers before implementation. Report generated on 2025-07-09 at 21:30 UTC.*

**CONFIDENTIAL - This document contains proprietary professional analysis.**