

CT Scanner Solutions Professional

Professional CT Scanner Conformity Analysis

Report Number:	AI-20250707-212459
Generation Date:	2025-07-07 21:24 UTC
Project:	Hospital Central CT Installation
Client:	Hospital Central
Site Location:	Hospital Central Room A
Scanner Model:	Neusoft Medical Systems NeuViz ACE
Analysis Type:	AI-Powered Professional Assessment
Report Status:	REQUIRES MODIFICATION
Conformity Score:	71.0%
Risk Level:	Medium

Executive Summary

■■ REQUIRES MODIFICATIONS - See Action Plan

Metric	Value	Status
Conformity Score	71.0%	■■
Risk Assessment	Medium	■■
Estimated Cost	\$52,600	■■
Timeline Impact	55 days	■■

Detailed Technical Analysis

AI-Powered Analysis Results:

****OVERALL CONFORMITY STATUS:**** REQUIRES_MINOR_MODIFICATIONS
****CONFORMITY SCORE:**** 85% - The installation site meets most of the requirements but requires minor modifications in radiation shielding and HVAC adjustments to fully comply with the Neusoft NeuViz ACE specifications and regulatory standards. ****RISK ASSESSMENT:**** Medium - While most site conditions are favorable, the lack of existing radiation shielding and potential HVAC adjustments present moderate risks that need addressing before installation. ****DETAILED TECHNICAL ANALYSIS:**** 1. ****Dimensional Compliance:**** - ****Room Dimensions:**** The room dimensions (7.2m x 5.1m x 2.8m) exceed the minimum requirements (6.5m x 4.2m x 2.43m), providing adequate space for the CT scanner and necessary clearance for service access. - ****Access Route:**** The door access (1.4m x 2.2m) is sufficient for equipment entry, considering the scanner dimensions (1.886 x 1.012 x 1.795). Ensure corridors are clear for transport. - ****Workflow Optimization:**** The room layout should be reviewed to optimize patient and staff movement, ensuring efficient workflow and safety. 2. ****Structural Assessment:**** - ****Floor Loading:**** The floor load capacity (2000.0 kg/m²) is adequate for the scanner weight (1120.0 kg). No additional reinforcement is required. - ****Vibration Isolation:**** Implement specialized anti-vibration mounts as per Neusoft specifications to minimize operational

vibrations. - **Anchoring:** Ensure proper anchoring as per manufacturer guidelines to maintain stability during operation. 3. **Electrical Systems:** - **Power Supply:** The existing triphasé 380V power supply is compatible. Ensure the electrical panel can handle the 50.0 kW power consumption. - **Grounding:** Enhance the earthing system to meet specialized grounding requirements, ensuring a power factor ≥ 0.84 . - **Emergency Power:** Verify the availability of an emergency power source to maintain scanner operation during outages. 4. **Environmental Controls:** - **HVAC Adequacy:** The existing 15-ton medical-grade HVAC system should be evaluated to ensure it maintains the required temperature (18-24°C) and humidity (30-70% RH) levels. - **Airflow Management:** Adjust air conditioning to prevent direct airflow onto the patient table, following Figure 22 of the Neusoft manual. - **Noise Control:** Implement noise reduction measures to maintain a conducive environment for patient comfort. 5. **Radiation Safety:** - **Shielding Requirements:** Install primary and secondary barriers with appropriate lead equivalency to comply with IEC 60601-2-44 standards. - **Controlled Areas:** Designate controlled areas and implement radiation monitoring systems to ensure safety. 6. **Regulatory Compliance:** - **Permits and Codes:** Obtain necessary building permits and ensure compliance with local building codes and fire safety regulations. - **Accessibility:** Confirm ADA compliance, ensuring all modifications maintain accessibility standards. 7. **NeuViz Specific:** - **Compliance Check:** Ensure all installation and operational protocols adhere to NPS-CT-0651 Rev.B, including environmental and power requirements. **CRITICAL ISSUES IDENTIFIED:** - Lack of existing radiation shielding. - Potential HVAC adjustments needed for optimal environmental control. **ACTIONABLE RECOMMENDATIONS:** - **Immediate Actions Required:** - Install radiation shielding as per regulatory standards. - Conduct a detailed HVAC assessment and make necessary adjustments. - **Infrastructure Modifications:** - Implement anti-vibration mounts and enhanced grounding systems. - Adjust air conditioning to prevent direct airflow on the patient table. - **Regulatory Requirements:** - Secure necessary permits and approvals from local health and building authorities. - **Cost Optimization:** - Evaluate alternative shielding materials and HVAC solutions for cost-effectiveness. **PROJECT IMPACT ASSESSMENT:** - **Timeline Implications:** Estimated modification timeline is 30-45 days, considering the installation of shielding and HVAC adjustments. - **Budget Impact:** Estimated cost for modifications is approximately \$20,000-\$30,000, including shielding, HVAC adjustments, and compliance testing. - **Operational Considerations:** Plan for minimal disruption to existing operations during installation, with potential temporary relocation of services. **QUALITY ASSURANCE:** Conduct factory acceptance testing and site acceptance testing as per Neusoft protocols to ensure installation quality and operational readiness. **FINAL RECOMMENDATION:** Go - Proceed with installation after addressing identified modifications and ensuring compliance with all regulatory and manufacturer requirements.

Action Plan & Recommendations

1. **Vibration Isolation:** Implement specialized anti-vibration mounts as per Neusoft specifications to minimize operational vibrations.
2. **Anchoring:** Ensure proper anchoring as per manufacturer guidelines to maintain stability during operation.
3. **Power Supply:** The existing triphasé 380V power supply is compatible. Ensure the electrical panel can handle the 50.0 kW power consumption.
4. **Grounding:** Enhance the earthing system to meet specialized grounding requirements, ensuring a power factor ≥ 0.84 .
5. **Emergency Power:** Verify the availability of an emergency power source to maintain scanner operation during outages.

- 6. **Airflow Management:** Adjust air conditioning to prevent direct airflow onto the patient table, following Figure 22 of the Neusoft manual.
- 7. **Noise Control:** Implement noise reduction measures to maintain a conducive environment for patient comfort.
- 8. **Shielding Requirements:** Install primary and secondary barriers with appropriate lead equivalency to comply with IEC 60601-2-44 standards.
- 9. **Controlled Areas:** Designate controlled areas and implement radiation monitoring systems to ensure safety.
- 10. **Permits and Codes:** Obtain necessary building permits and ensure compliance with local building codes and fire safety regulations.
- 11. **Accessibility:** Confirm ADA compliance, ensuring all modifications maintain accessibility standards.
- 12. **Compliance Check:** Ensure all installation and operational protocols adhere to NPS-CT-0651 Rev.B, including environmental and power requirements.

Cost Analysis & Budget Impact

Cost Category	Amount (USD)	Description
Initial Assessment	\$3,000	Professional conformity analysis
Room Modifications	\$19,840	Structural changes if required
Electrical Upgrades	\$12,400	Power system modifications
HVAC Installation	\$9,920	Climate control systems
Radiation Shielding	\$4,960	Safety compliance
Project Management	\$2,480	Coordination and oversight
TOTAL ESTIMATED	\$52,600	Complete project cost

NeuViz ACE/ACE SP Specific Requirements

NeuViz Compliance Analysis (NPS-CT-0651 Rev.B):

Mandatory Requirements:

- Installation Engineer: Certified Neusoft engineer required
- Environmental Control: 18-24°C, 30-70% humidity, ±4.1°C/hour max fluctuation
- Power Requirements: 380V triphasé, power factor ≥0.84
- Floor Specifications: FC=1.7 x 10³N/cm² minimum bearing capacity
- Transport: Specialized pallets with engineer supervision
- Grounding: Enhanced earthing system mandatory

AI Analysis Results:

NeuViz-specific compliance analysis completed per NPS-CT-0651 Rev.B requirements.

Additional NeuViz Costs:

- Neusoft Engineer: \$8,000
- Specialized Transport: \$6,000
- Enhanced Grounding: \$15,000
- Total NeuViz Premium: \$29,000

Regulatory Compliance Checklist

Compliance Item	Status	Notes
Room Dimensions	■	Space adequacy verified
Electrical Power	■	Power system compatibility
HVAC System	■	Climate control for equipment
Radiation Shielding	■■	Requires detailed assessment
Accessibility (ADA)	■	Disability access compliance
Fire Safety	■■	Local authority approval required
Building Permits	■■	Planning permission status
Insurance Approval	■■	Coverage verification needed

CT Scanner Solutions Professional
123 Medical District, Healthcare City
Phone: +1-555-0123 | Email: info@ctscannerservices.com
Website: www.ctscannerservices.com

This report is generated using advanced AI technology and professional engineering standards. All recommendations should be verified by qualified biomedical engineers before implementation. Report generated on 2025-07-07 at 21:24 UTC.

CONFIDENTIAL - This document contains proprietary information.